



101 Tosca Drive  
Stoughton, MA 02072 USA  
(Phone) 781.297-2034  
(FAX) 781.297-2039  
(Web) [www.csmisolutions.com](http://www.csmisolutions.com)

Congratulations on your decision to select the **HUMAC® NORM™ Testing & Rehabilitation System** for your facility. We appreciate your business and promise to do our best to make sure that your order, shipment, delivery and installation are handled to the best of our ability.

One of the key elements to our most successful installations is facility planning. This Site Survey has been developed over the years to assist in completing your order with minimum inconvenience and business interruptions. It is important that the completed survey be returned in order for CSMi to process your NORM order. Upon completion, please email [info@csmisolutions.com](mailto:info@csmisolutions.com) or FAX 781-297-2039 the completed copy to CSMi. If you have any questions, call Rob at 781-818-3804.

Thank you in advance for your assistance.

**IMPORTANT NOTE:**

- 1. Please make sure your equipment is delivered to the exact location of its use and the location and access match the answers below. There will be an additional charge if CSMi is required to re-visit the facility or move your equipment to complete the installation*
- 2. A receiving dock is not required as our trucks have a lift-gate.*

<b>FACILITY NAME:</b>	
<b>DEPARTMENT:</b>	
<b>ADDRESS:</b>	
<b>CITY, STATE, ZIP:</b>	
<b>COUNTRY:</b>	
<b>PRIMARY CONTACT:</b>	
<b>PC'S PHONE:</b>	
<b>PC'S EMAIL:</b>	
<b>ALTERNATE CONTACT:</b>	
<b>AC'S PHONE:</b>	
<b>AC'S EMAIL:</b>	
<b>FORM COMPLETED BY:</b>	
<b>SIGNATURE:</b>	

## FACILITY

1. Is your facility/dept. complete and ready to accept delivery of your NORM? ☐ Yes ☐ No

If not, when would you be ready to accept a delivery? Date \_\_\_\_\_

2. Has your flooring and /or carpet been installed? ☐ Yes ☐ No

What type of flooring/carpeting do you have? \_\_\_\_\_

3. What are the dimensions of the room where the NORM will be located?

Height: \_\_\_\_\_ Width: \_\_\_\_\_ Length: \_\_\_\_\_

4. On which floor of the building will the NORM be placed? \_\_\_\_\_

Is there an elevator? ☐ Yes ☐ No Width of elevator door (inches)? \_\_\_\_\_

Are there stairs? ☐ Yes ☐ No How many flights? \_\_\_\_\_

Are there any unusual or narrow corner ways or entrances? ☐ Yes ☐ No

Please explain: \_\_\_\_\_

What is the width of narrowest doorway through which the NORM will be moved? \_\_\_\_\_

5. Do you have a receiving dock? ☐ Yes ☐ No

6. Is it standard Tailboard height? ☐ Yes ☐ No

If not, how do you receive equipment? \_\_\_\_\_

7. How many feet of operating space have been allocated for the use of the NORM (Refer to FLOOR SPACE REQUIREMENTS on the last page)? \_\_\_\_\_X\_\_\_\_\_

8. Will an internet connection to the NORM computer be available? ☐ Yes ☐ No

**Note:** This is recommended but not required. An internet connection allows CSMi to provide enhanced remote product support.

## ELECTRICAL REQUIREMENTS

1. Have you installed an independent, dedicated 208VAC (Nominal), 50/60 Hz, Single Phase, 20 AMP line, no more than 10 feet from the Isolation Transformer? ☐ Yes ☐ No

The Isolation Transformer is provided by CSMi and included with HUMAC NORM.

**NOTE:** The HUMAC/NORM System is compatible with:

200/208/220/230/240VAC ~,

50/60 HZ

Single-Phase 1Ø.


2. Is the HUMAC NORM's electrical receptacle a NEMA number 6-20R Isolated Ground (Hubbell IG5461) or equivalent Hospital Grade (Hubbell 8410IG); does it match the diagram below? ☐ Yes ☐ No



<b>Rating:</b>	250Volts, 20 Amps.
<b>Contact:</b>	Solid Brass.
<b>Installation:</b>	The receptacle should be mounted with the <b>ground positioned at the top.</b>
<b>Color:</b>	Orange is recommended but not required. Refer to local building ordinances for specific requirements.

3. Is all electrical work in compliance with the local building codes and all other agencies having jurisdiction? ☐ Yes ☐ No

## SUGGESTED FLOOR PLAN



063-1

Figure 1. Suggested Floor Plan.